KONICA MINOLTA

REGIUS CASSETTE/PLATE

CAUTION Cautions Related to Cassette Handling

EN

- When an artifact is found on the exposed image check the plate. If any scratches, dust or foreign substances are present on the plate, clean plate. After cleaning, perform a test exposure to check if problem has been resolved. If the artifact is still present causing difficulty in medical diagnosis, stop using the cassette and replace with new.
- To prevent cross contamination always clean the cassette according to proper cleaning instructions.
- If patient's blood and/or bodily fluids are present on the cassette, repair may be refused. Always clean the cassette appropriately prior to sending them in for repair.
- Dropping, bending or applying a shock to the cassette may lead to deformation or breakage. Always handle them with care. If the cassette has been dropped, bent, or has been subjected to excessive force, verify the cassette has no chips, cracks, or other defects. Make sure to carefully examine the corners and the lock of the cassette before use. Use of faulty cassettes can cause injury or cassette transport error. If damage is observed, immediately stop using the cassette.
- Read precautions detailed in the operation manual of the Reader prior to handling cassettes.

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Konica Minolta Healthcare Americas, Inc.
U.S.A.
Konica Minolta Business Solutions Europe GmbH
Konica Minolta Business Solutions (Canada) Ltd.
Konica Minolta Medical & Graphic (SHANGHAI) Co., Ltd.
Konica Minolta Healthcare India Private Ltd.