



KONICA MINOLTA

MEDICAL IMAGING FILM

SR-IC PLUS

I · M · A · G · I · N · G




Giving Shape to Ideas

MEDICAL IMAGING FILM SR-IC PLUS

*Produces images of high-quality and sharp detail over a wide density range.
Will not miss even minute details for diagnosis.*

Characteristic curve



Description

KONICAMINOLTA SR-IC PLUS is an orthochromatic medical imaging film designed with excellent contrast, gradation and sharp detail with minimal base fog due to the fine grain emulsion on one side and an anti-halation layer on the other.

Application

M.R.I./CT scanning/D.S.A./Ultra sound/Nuclear imaging/
Thermography/Mammography

*Subject to change without further notice.

Konica Minolta medical film is manufactured with a quality management system that has been certified to be in compliance with the ISO 9001:2000 and ISO 13485:2003 quality management standards, as well as with the medical device directive (MDD) 93/42/EEC.



KONICA MINOLTA

KONICA MINOLTA, INC.

1 Sakura-machi, Hino-shi, Tokyo, 191-8511, Japan

Distributed by :