



KONICA MINOLTA

PRIMARY IMAGING SOLUTIONS

Konica Minolta. Right Solutions. Right Time.

REGIUS Xpress & Nano CR

Superior image quality, fast and reliable, easy to use

The REGIUS Xpress and Nano Computed Radiography systems offer superior image quality, are fast and reliable, and easy to use. They have earned the highest customer satisfaction rating in the industry.

- Superior image quality is a hallmark of Konica Minolta CR systems. All Xpress and Nano models include patented hybrid processing that divides an image into multiple frequency bands producing a more natural looking image with less noise and artifacts.
- Fast and reliable imaging solutions are imperative in today's health care environment. Xpress CR offers the high performance of up to 92 plates per hour, while Nano CR can process up to 76 plates per hour for fast patient turnaround and improved workflow. Our reliability is the highest customer satisfaction rating in the industry.
- Konica Minolta CR systems are simple and easy to use. Our innovative technology is designed to make your job easier. Touch screen control panels provide intuitive software tools for generating superior images in a few easy steps.



Superior Image Quality

Clinicians need high-quality images to make clinically confident patient care decisions. Konica Minolta is a leader in delivering the highest quality images across every imaging application. From the private practice to the enterprise, you can rely on Konica Minolta's patented, hybrid processing to deliver superior images with the clinical information you need in each patient exposure.

Giving Shape to Ideas

REGIUS Xpress & Nano CR

*Superior image quality, fast
and reliable, easy to use*



REGIUS Model 210

Multi Bay CR Reader

Models:

Xpress CR 'E' Class

Xpress CR 'A' Class

For Use In:

Enterprise Imaging

Busy Imaging Centers

Large Orthopedic Centers

NANO CR

REGIUS Model 110

Single Bay CR Reader

Models:

Nano CR 'E' Class

Nano CR 'A' Class

Nano CR 'P+' Class (Private Practice)

For Use In:

Community Hospitals

Imaging Centers

Orthopedic Centers

Fast and Reliable

With a fast, efficient workflow and the reliability to consistently deliver superior imaging across exams and patient types, the Konica Minolta CR products offer high economic value in a cost-constrained healthcare environment. Xpress CR generates 92 plates per hour while the Nano CR offers up to 76 plates per hour so your department can meet high to mid-volume needs.

Easy to use

You need to stay focused on your patients, not the equipment. You can depend on Konica Minolta for CR imaging that is intelligent and easy-to-use, even for the novice technologist. Intelligent options for full user functionality, DICOM store and print options, and quick image processing provide worry-free operation every time.

REGIUS CR Systems are Preferred Products for:

Enterprise Hospital, Imaging Centers, Orthopedic Centers, Small & Large Clinics, Veterinary, Chiropractic and Private Practice facilities.

REGIUS CR Systems Offer:

Factory Refresh packages to upgrade CS-1/CS-2/CS-3 (XP based) Control Stations to CS-7 (Windows 7) Control Stations to help meet today's HIPAA and/or facility requirements.



KONICA MINOLTA

Konica Minolta Medical Imaging

411 Newark Pompton Turnpike

Wayne, NJ 07470

Tel: (973) 633-1500 Fax: (973) 523-7408

konicaminolta.com/medicalusa

REGIUS Xpress & Nano CR 2/2015