




***Produces images of high-quality and sharp detail over a wide density range. Reproduces CRT images with high contrast and resolution.***

**Characteristic Curve****Description**

KONICA MINOLTA SR-IC PLUS is an orthochromatic video imaging film designed with excellent contrast, gradation and sharp detail with minimal base fog due to the fine grain emulsion on one side and an anti-halation layer on the other. It may be used in a super rapid processing cycle.

**Application**

M.R.I./CT scanning/D.S.A./Ultra sound/Nuclear imaging/Thermography/Mammography

\*Subject to change without further notice.

Konica Minolta medical film is manufactured with a quality management system that has been certified to be in compliance with the ISO 9001:2000 and ISO 13485:2003 quality management standards, as well as with the medical device directive (MDD) 93/42/EEC.

**KONICA MINOLTA****KONICA MINOLTA MEDICAL IMAGING USA, INC.**

411 NEWARK POMPTON TURNPIKE, WAYNE, NJ, 07470, USA

Tel: (973) 633-1500 Fax: (973) 633-0562, [medical.konicaminolta.us](http://medical.konicaminolta.us)

Distributed by: